

beta galactosidase air ih sample method

shipping requirements	On ice pack
sampling media	25 mm Teflon Filter, 1 micron pore size (Millipore Cat. # FALP02500); cellulose support pad, (SKC Cat. # 225-28); in 2-pc polypropylene cassette, (SKC Cat# 225-8585)
sampling rate	2 L/min
recommended minimum sampling volume	50 L

If you wish to submit a sample, please follow the sample submission procedures found at safebridge.com/laboratory/sample-submission. SafeBridge regularly develops new and custom industrial hygiene air and surface sampling methods. For more information, please call (650) 961-4820. SafeBridge Consultants, Inc. is the premier occupational health and safety and risk assessment consulting firm to the pharmaceutical, biotechnology, and fine chemical industry. Learn more at safebridge.com.