

chondroitinase surface ih sample method

shipping requirements	On ice pack, overnight
swabbing media	Two CleanTips swabs TX714A; in polypropylene test tube with dual position cap containing 1 mL extracting buffer
swabbing solution	0.05% tween 20 in 0.05 M carbonate/bicarbonate buffer, pH 9.6

If you wish to submit a sample, please follow the sample submission procedures found at safebridge.com/laboratory/sample-submission. SafeBridge regularly develops new and custom industrial hygiene air and surface sampling methods. For more information, please call (650) 961-4820. SafeBridge Consultants, Inc. is the premier occupational health and safety and risk assessment consulting firm to the pharmaceutical, biotechnology, and fine chemical industry. Learn more at safebridge.com.