

insulin air ih sample method

storage conditions	At least 2 weeks at ~4 °C
shipping requirements	Cold, dry ice
recommended minimum sampling volume	120L

If you wish to submit a sample, please follow the sample submission procedures found at safebridge.com/laboratory/sample-submission.

SafeBridge regularly develops new and custom industrial hygiene air and surface sampling methods. For more information, please call (650) 961-4820.

SafeBridge Consultants, Inc. is the premier occupational health and safety and risk assessment consulting firm to the pharmaceutical, biotechnology, and fine chemical industry. Learn more at safebridge.com.