

insulin porcine surface ih sample method

shipping requirements	Ambient temperature, overnight
swabbing media	Two CleanTips Swabs, TX740B with CleanFoam (polyurethane) head; in polypropylene test tube with dual position cap
swabbing solution	5% Acetic acid/DI Water

If you wish to submit a sample, please follow the sample submission procedures found at safebridge.com/laboratory/sample-submission. SafeBridge regularly develops new and custom industrial hygiene air and surface sampling methods. For more information, please call (650) 961-4820. SafeBridge Consultants, Inc. is the premier occupational health and safety and risk assessment consulting firm to the pharmaceutical, biotechnology, and fine chemical industry. Learn more at safebridge.com.