

irinotecan hcl surface ih sample method

storage conditions	At ~ ?20 ? for at least 21 days, protected from light
shipping requirements	Ice pack, protected from light
swabbing media	Two CleanTips swabs, TX714A
swabbing solution	0.1% TFA in DI water

If you wish to submit a sample, please follow the sample submission procedures found at safebridge.com/laboratory/sample-submission. SafeBridge regularly develops new and custom industrial hygiene air and surface sampling methods. For more information, please call (650) 961-4820. SafeBridge Consultants, Inc. is the premier occupational health and safety and risk assessment consulting firm to the pharmaceutical, biotechnology, and fine chemical industry. Learn more at safebridge.com.