

spinetoram surface ih sample method

shipping requirements	Ambient temperature
swabbing media	3 Large Alpha swabs, CleanTip swab, cat. TX714A, TWTexwipe
swabbing solution	Methanol: DI water (80:20, v/v)

If you wish to submit a sample, please follow the sample submission procedures found at safebridge.com/laboratory/sample-submission. SafeBridge regularly develops new and custom industrial hygiene air and surface sampling methods. For more information, please call (650) 961-4820. SafeBridge Consultants, Inc. is the premier occupational health and safety and risk assessment consulting firm to the pharmaceutical, biotechnology, and fine chemical industry. Learn more at safebridge.com.