

# Tretinoin Air IH Sample Method

|                                            |                                                                                       |
|--------------------------------------------|---------------------------------------------------------------------------------------|
| <b>shipping requirements</b>               | Cold on ice pack, protected from light                                                |
| <b>sampling media</b>                      | 37 mm teflon filter, 1 micron pore size, support pad, in 2-piece polystyrene cassette |
| <b>sampling rate</b>                       | 3 L/min                                                                               |
| <b>recommended minimum sampling volume</b> | 100 L                                                                                 |

Generated on October 1, 2025